



## Single Fold STANAG 2040 NATO Stretcher

Single Fold STANAG 2040 Stretcher by Golden Season Pte Ltd



Golden Season's Single Fold STANAG 2040 Stretcher serves as a reliable casualty evacuation platform and is manufactured in accordance with NATO Standards (STANAG 2040).



This appliance conforms to the Directive 93/42/CEE "Medical Devices"

Technical Specifications	
Dimensions Extended	2290mm X 570mm X 153mm
Dimensions Folded (Handles Retracted)	1980mm X 85mm X 153mm
Canvas Specification	Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards
Frame	Aerospace Grade Aluminum
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7.5kg
CE Mark	Yes

### Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season Building) T 65.6863 2228 F 65.6363 2478  
Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) Co. Regn No: 03734/1984E  
GST Registration No. M2-0067495-2