



Roller Stair Chair



The Roller Stair Chair is designed for transporting a patient from multi-level buildings safely and efficiently. Equipped with a set of rollers, the Roller Stair Chair enables a single user to “glide” down a flight of stairs. The Roller Stair Chair is foldable for compact storage, and easily deployable. The Roller Stair Chair can also function as a normal wheelchair.

Specifications	
Unfolded size:	98 x 54 x 80cm
Folded size:	98 x 54 x 35cm
Net weight:	15kg
Load Bearing:	200kg

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd