



Scoop Stretcher



The Golden Season Scoop Stretcher is designed for transportation of patients suffering from back injuries without adding additional trauma. Standard evacuation platforms contact the patient's spine adding pain and discomfort during transportation.

Scooping up the patient

The scoop stretcher is able to effectively scoop up the patient from laying position.

Foldable design

Our scoop stretchers are foldable and compact for storage. It is easy to transport as well.

Extremely durable

Golden Season's scoop the stretchers utilize rugged joints and an adamant construction to with stand abuse and rough handling.

Specifications:	
Dimension (Deploy):	210cm x 44cm x 9cm
Dimension (Folded):	120cm x 44cm x 9cm
Weight:	Approx. 7.5kg
Accessories:	2 / 3 patient restraints provided

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd