



Stair Chair Deluxe



The Stair Chair Deluxe assist paramedics and health care professionals to transport the patient up or down the stairways. It's also comes with equipped with 4 wheels instead of 2.

Compactible Design

The Stair Chair Deluxe is compactible for easy and deploys in an instant.

Ergonomic retractable Handles

Golden Season's Stair Chair Deluxe is designed with ergonomic retractable handles for the comfort of health care professionals in executing their duties.

Come with Patient Restraints

Stair Chair Deluxe comes with patient restraints to strap the patient to the chair , enhancing the safety of the patient during transport.

Specifications	
Dimension when Extended:	120cm x 50cm x 93.5cm
Dimension when Folded:	23cm x 50cm x 93.5cm
Weight:	10.13kg
Number of Wheels:	4 Wheels
Accessories:	Patient Restraints (2 sets are provided)

*Specifications are subject to change without prior notice.

Golden Season Pte Ltd