

## **Stretcher Wheel Cart/Wheeled Litter Carrier**

(PATENTED)



Golden Season's Stretcher Wheel Cart/Wheeled Litter Carrier is a lightweight casualty transport system that allows one person to transport a casualty strapped to a stretcher, and also performs as a litter stand in the battlefield or medical stations.

The revolutionary design is foldable, and can also be totally dissembled for compact storage, easy shipping and transportation. The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern combat litters, and is equipped with non puncture wheels. A single user can operate the Stretcher Wheel Cart/Wheeled Litter Carrier.



| Technical Information   |                       |
|-------------------------|-----------------------|
| Deployed dimensions:    | L143 x W70 x H92 (cm) |
| Half Folded dimensions: | L126 x W70 x H50 (cm) |

## **Golden Season Pte Ltd**

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478 Email: info@goldenseason.com.sg | www.gseason.com Co. Regn No: 03734/1984E GST Registration No. M2-0067495-2

## **GOLDEN SEASON** PTE. LTD.

www.goldenseason.com.sg | www.gseason.com

| Fully Folded dimensions: | L107 x W70 x H25 (cm)    |
|--------------------------|--------------------------|
| Weight:                  | 17.5kgs excluding litter |

\*Specifications are subject to change without prior notice.



## **Characteristics**

- ✓ Foldable and compactable
- $\checkmark$  Able to perform as a litter stand
- ✓ Lightweight 17.5kgs only
- ✓ Pneumatic non puncture wheels
- ✓ Can be fully disassembled for easy shipping

The Stretcher Wheel Cart/Wheeled Litter Carrier is compatible with modern stretcher systems including the Medevac4, TalonII 90C, Raven 90C, T3 Assault Litters and other stretchers.

**Golden Season Pte Ltd**