



## Emergency Trauma Dressing (ETD) - 4 in.



The Emergency Trauma Dressing (ETD™) series of bandages are engineered using data from over 63 After Action Reviews (AARs) of real-world battle dressing applications. The ETD™ consists of a resilient elastic wrap equipped with a sterile non-adherent pad, durable securing device and innovative Quick-Grip™ Roll Control configuration. This combination facilitates a rapid, controlled, and extremely effective application process in the most demanding conditions. The ETD™ can also be used to sling and swath an arm, to secure splints to an extremity or to perform in any other function requiring an elastic wrap.

### Special Features:

- ✓ Sterile bandage with up to 5 year expiry, Made in the U.S. & Berry Amendment Compliant
- ✓ Suitable for applying immediate direct pressure
- ✓ Multi-functional, all-in-one design
- ✓ Simple application that is easy to remember under combat/survival stress
- ✓ Omission of pressure bars and hooks to moderate the use of fine motor skills
- ✓ Durable securing device
- ✓ Quick-Grip™ Roll Control tabs to prevent the bandage from unrolling during application
- ✓ Vacuum-sealed, Low-Cube Packaging
- ✓ Cost effective

# GOLDEN SEASON PTE. LTD.

[www.goldenseason.com.sg](http://www.goldenseason.com.sg) | [www.gseason.com](http://www.gseason.com)



## Contents:

Item#: 30-0031 NSN#: 6510-01-558-4108

- ✓ Sterile bandage with up to 5 year expiry, Made in the U.S. & Berry Amendment Compliant
- ✓ Suitable for applying immediate direct pressure
- ✓ Simple application that is easy to remember under combat/survival stress
- ✓ Omission of pressure bars and hooks to moderate the use of fine motor skills
- ✓ Quick-Grip™ Roll Control tabs to prevent the bandage from unrolling during application
- ✓ Vacuum-sealed, Low-Cube Packaging

Specifications	
Packaged (4 in.):	L 4.5 in. x W 2.5 in. x D 1.75 in.
Bandage (4 in.):	L 70 in. x W 4 in.
Weight (4 in.):	2.6 oz

\*Specifications are subject to change without prior notice.

---

## Golden Season Pte Ltd

31 Playfair Road Singapore 367993 (Golden Season building) T 65.6863 2228 F 65.6363 2478  
Email: [info@goldenseason.com.sg](mailto:info@goldenseason.com.sg) Website: [www.goldenseason.com.sg](http://www.goldenseason.com.sg) | [www.gseason.com](http://www.gseason.com) Co. Regn No: 03734/1984E  
GST Registration No. M2-0067495-2