



## Single Fold STANAG 2040 NATO Stretcher



Golden Season's Single Fold STANAG 2040 Stretcher serves as a reliable casualty evacuation platform and is manufactured in accordance with NATO Standards (STANAG 2040).

**CE** This appliance conforms to the Directive 93/42/CEE "Medical Devices"

Specifications:	
Dimensions Extended	2290mm X 570mm X 153mm
Dimensions Folded (Handles Retracted)	1980mm X 85mm X 153mm
Canvas Specification	Warp: Minimum 2400N / 50mm Weft: Minimum 1400N / 50mm Compliant with MIL-L-16446A Standards
Frame	Aerospace Grade Aluminum
Handles	Retractable, Automatic locking, Ergonomically Designed
Weight	7.5kg
CE Mark	Yes

\*Specifications are subject to change without prior notice.

### Golden Season Pte Ltd